

Anti-Phospho-CDCP1-Tyr806 antibody (STJ196272)

STJ196272

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Cub Domain-Containing Protein 1-Tyr806 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64866
Gene Symbol	CDCP1
Uniprot ID	CDCP1_HUMAN
Immunogen	Synthesized phosho peptide around human CDCP1 (Tyr806)
Region	
Specificity	Phospho-CDCP1-Tyr806 polyclonal antibody (Cub Domain-Containing Protein 1) binds to endogenous Cub Domain-Containing Protein 1 only when phosphorylated at Tyr806.
Immunogen Sequence	