

Anti-Phospho-KIT-Tyr568/570 antibody (STJ196274)

STJ196274

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Mast/Stem Cell Growth Factor Receptor Kit-Tyr568/570 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3815
Gene Symbol	KIT
Uniprot ID	KIT_HUMAN
Immunogen	Synthesized phosho peptide around human c-Kit (Tyr568 and 570)
Region	
Specificity	Phospho-KIT-Tyr568/570 polyclonal antibody (Mast/Stem Cell Growth Factor Receptor Kit) binds to endogenous Mast/Stem Cell Growth Factor Receptor Kit only when phosphorylated at Tyr568/570.
Immunogen Sequence	