

Anti-Phospho-CTDSPL2-Ser104 antibody (STJ196277)

STJ196277

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Ctd Small Phosphatase-Like Protein 2-Ser104 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51496
Gene Symbol	CTDSPL2
Uniprot ID	CTSL2_HUMAN
Immunogen	Synthesized phosho peptide around human CTDSPL2 (Ser104)
Immunogen Region	
Specificity	Phospho-CTDSPL2-Ser104 polyclonal antibody (Ctd Small Phosphatase-Like Protein 2) binds to endogenous Ctd Small Phosphatase-Like Protein 2 only when phosphorylated at Ser104.
Immunogen Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081