

Anti-Phospho-CCAR2-Thr454 antibody (STJ196283)

STJ196283

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Cell Cycle And Apoptosis Regulator Protein 2-Thr454 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57805
Gene Symbol	CCAR2
Uniprot ID	CCAR2_HUMAN
Immunogen	Synthesized phosho peptide around human DBC1 (Thr454)
Region	
Specificity	Phospho-CCAR2-Thr454 polyclonal antibody (Cell Cycle And Apoptosis Regulator Protein 2) binds to endogenous Cell Cycle And Apoptosis Regulator Protein 2 only when phosphorylated at Thr454.
Immunogen Sequence	