

Anti-Phospho-DSP-Ser165/166 antibody (STJ196285)

STJ196285

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Desmoplakin-Ser165/166 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1832
Gene Symbol	DSP
Uniprot ID	DESP_HUMAN
Immunogen	Synthesized phospho peptide around human Desmoplakin (Ser165 and 166)
Immunogen Region	
Specificity	Phospho-DSP-Ser165/166 polyclonal antibody (Desmoplakin) binds to endogenous Desmoplakin only when phosphorylated at Ser165/166.
Immunogen Sequence	