

Anti-Phospho-FOXO3-Ser294 antibody (STJ196308)

STJ196308

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Forkhead Box Protein O3-Ser294 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2309
Gene Symbol	FOXO3
Uniprot ID	FOXO3_HUMAN
Immunogen	Synthesized phosho peptide around human FoxO3a (Ser294)
Immunogen Region	
Specificity	Phospho-FOXO3-Ser294 polyclonal antibody (Forkhead Box Protein O3) binds to endogenous Forkhead Box Protein O3 only when phosphorylated at Ser294.
Immunogen Sequence	