

Anti-Phospho-GRB10-Ser476 antibody (STJ196318)

STJ196318

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Growth Factor Receptor-Bound Protein 10-Ser476 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2887
Gene Symbol	GRB10
Uniprot ID	GRB10_HUMAN
Immunogen	Synthesized phospho peptide around human GRB10 (Ser476)
Immunogen Region	
Specificity	Phospho-GRB10-Ser476 polyclonal antibody (Growth Factor Receptor-Bound Protein 10) binds to endogenous Growth Factor Receptor-Bound Protein 10 only when phosphorylated at Ser476.
Immunogen Sequence	