

Anti-Phospho-GSK3B-Thr390 antibody (STJ196319)

STJ196319

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Glycogen Synthase Kinase-3 Beta-Thr390 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2932
Gene Symbol	GSK3B
Uniprot ID	GSK3B_HUMAN
Immunogen	Synthesized phosho peptide around human GSK-3 Beta (Thr390)
Region	
Specificity	Phospho-GSK3B-Thr390 polyclonal antibody (Glycogen Synthase Kinase-3 Beta) binds to endogenous Glycogen Synthase Kinase-3 Beta only when phosphorylated at Thr390.
Immunogen Sequence	