

Anti-Phospho-ERBB4-Tyr984 antibody (STJ196321)

STJ196321

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Receptor Tyrosine-Protein Kinase Erbb-4-Tyr984 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2066
Gene Symbol	ERBB4
Uniprot ID	ERBB4_HUMAN
Immunogen	Synthesized phospho peptide around human HER4 and ErbB4 (Tyr984)
Immunogen Region	
Specificity	Phospho-ERBB4-Tyr984 polyclonal antibody (Receptor Tyrosine-Protein Kinase Erbb-4) binds to endogenous Receptor Tyrosine-Protein Kinase Erbb-4 only when phosphorylated at Tyr984.
Immunogen Sequence	