

Anti-Phospho-HSP90AA1-Thr5/7 antibody (STJ196324)

STJ196324

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Heat Shock Protein Hsp 90-Alpha-Thr5/7 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3320
Gene Symbol	HSP90AA1
Uniprot ID	HS90A_HUMAN
Immunogen	Synthesized phospho peptide around human HSP90 Alpha (Thr5 and 7)
Immunogen Region	
Specificity	Phospho-HSP90AA1-Thr5/7 polyclonal antibody (Heat Shock Protein Hsp 90-Alpha) binds to endogenous Heat Shock Protein Hsp 90-Alpha only when phosphorylated at Thr5/7.
Immunogen Sequence	