

Anti-Phospho-IGF1R-Tyr1131 antibody (STJ196325)

STJ196325

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Insulin-Like Growth Factor 1 Receptor-Tyr1131 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3480
Gene Symbol	IGF1R
Uniprot ID	IGF1R_HUMAN
Immunogen	Synthesized phosho peptide around human IGF-I Receptor Beta (Tyr1131)
Immunogen Region	
Specificity	Phospho-IGF1R-Tyr1131 polyclonal antibody (Insulin-Like Growth Factor 1 Receptor) binds to endogenous Insulin-Like Growth Factor 1 Receptor only when phosphorylated at Tyr1131.
Immunogen Sequence	