

Anti-Phospho-IRF7-Ser477 antibody (STJ196333)

STJ196333

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Interferon Regulatory Factor 7-Ser477 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3665
Gene Symbol	IRF7
Uniprot ID	IRF7_HUMAN
Immunogen	Synthesized phospho peptide around human IRF-7 (Ser477)
Immunogen Region	
Specificity	Phospho-IRF7-Ser477 polyclonal antibody (Interferon Regulatory Factor 7) binds to endogenous Interferon Regulatory Factor 7 only when phosphorylated at Ser477.
Immunogen Sequence	