

Anti-Phospho-JAK1-Tyr1034/1035 antibody (STJ196339)

STJ196339

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Tyrosine-Protein Kinase Jak1-Tyr1034/1035 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3716
Gene Symbol	JAK1
Uniprot ID	JAK1_HUMAN
Immunogen	Synthesized phosho peptide around human Jak1 (Tyr1034 and 1035)
Immunogen Region	
Specificity	Phospho-JAK1-Tyr1034/1035 polyclonal antibody (Tyrosine-Protein Kinase Jak1) binds to endogenous Tyrosine-Protein Kinase Jak1 only when phosphorylated at Tyr1034/1035.
Immunogen Sequence	