

Anti-Phospho-JAK2-Tyr1008 antibody (STJ196340)

STJ196340

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Tyrosine-Protein Kinase Jak2-Tyr1008 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3717
Gene Symbol	JAK2
Uniprot ID	JAK2_HUMAN
Immunogen	Synthesized phosho peptide around human Jak2 (Tyr1008)
Immunogen Region	
Specificity	Phospho-JAK2-Tyr1008 polyclonal antibody (Tyrosine-Protein Kinase Jak2) binds to endogenous Tyrosine-Protein Kinase Jak2 only when phosphorylated at Tyr1008.
Immunogen Sequence	