

Anti-Phospho-JAK3-Tyr980/981 antibody (STJ196341)

STJ196341

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Tyrosine-Protein Kinase Jak3-Tyr980/981 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3718
Gene Symbol	JAK3
Uniprot ID	JAK3_HUMAN
Immunogen	Synthesized phosho peptide around human Jak3 (Tyr980 and 981)
Immunogen Region	
Specificity	Phospho-JAK3-Tyr980/981 polyclonal antibody (Tyrosine-Protein Kinase Jak3) binds to endogenous Tyrosine-Protein Kinase Jak3 only when phosphorylated at Tyr980/981.
Immunogen Sequence	