

Anti-Phospho-LYN-Tyr507 antibody (STJ196350)

STJ196350

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Tyrosine-Protein Kinase Lyn-Tyr507 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4067
Gene Symbol	LYN
Uniprot ID	LYN_HUMAN
Immunogen	Synthesized phosho peptide around human Lyn (Tyr507)
Immunogen Region	
Specificity	Phospho-LYN-Tyr507 polyclonal antibody (Tyrosine-Protein Kinase Lyn) binds to endogenous Tyrosine-Protein Kinase Lyn only when phosphorylated at Tyr507.
Immunogen Sequence	