

Anti-Phospho-NEDD4L-Ser342 antibody (STJ196363)

STJ196363

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-E3 Ubiquitin-Protein Ligase Nedd4-Like-Ser342 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23327
Gene Symbol	NEDD4L
Uniprot ID	NED4L_HUMAN
Immunogen	Synthesized phospho peptide around human NEDD4L (Ser342)
Immunogen Region	
Specificity	Phospho-NEDD4L-Ser342 polyclonal antibody (E3 Ubiquitin-Protein Ligase Nedd4-Like) binds to endogenous E3 Ubiquitin-Protein Ligase Nedd4-Like only when phosphorylated at Ser342.
Immunogen Sequence	