

Anti-Phospho-NEDD4L-Ser448 antibody (STJ196364)

STJ196364

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-Phospho-E3 Ubiquitin-Protein Ligase Nedd4-Like-Ser448 is suitable for use in Western Blot research applications. |
| Applications | WB |
| Host/Source | Rabbit |
| Reactivity | Human, Rat, Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography. |
| Dilution Range | WB 1:1000-2000 |
| Formulation | PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide. |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 23327 |
| Gene Symbol | NEDD4L |
| Uniprot ID | NED4L_HUMAN |
| Immunogen | Synthesized phospho peptide around human NEDD4L (Ser448) |
| Immunogen Region | |
| Specificity | Phospho-NEDD4L-Ser448 polyclonal antibody (E3 Ubiquitin-Protein Ligase Nedd4-Like) binds to endogenous E3 Ubiquitin-Protein Ligase Nedd4-Like only when phosphorylated at Ser448. |
| Immunogen Sequence | |