

Anti-Phospho-PAK1-Ser20 antibody (STJ196377)

STJ196377

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Serine/Threonine-Protein Kinase Pak 1-Ser20 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5058
Gene Symbol	PAK1
Uniprot ID	PAK1_HUMAN
Immunogen	Synthesized phospho peptide around human PAK2 (Ser20)
Immunogen Region	
Specificity	Phospho-PAK1-Ser20 polyclonal antibody (Serine/Threonine-Protein Kinase Pak 1) binds to endogenous Serine/Threonine-Protein Kinase Pak 1 only when phosphorylated at Ser20.
Immunogen Sequence	