

Anti-Phospho-PDGFR α -Tyr1018 antibody (STJ196378)

STJ196378

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Platelet-Derived Growth Factor Receptor Alpha-Tyr1018 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5156
Gene Symbol	PDGFRA
Uniprot ID	PGFRA_HUMAN
Immunogen	Synthesized phosho peptide around human PDGF Receptor Alpha (Tyr1018)
Region	
Specificity	Phospho-PDGFR α -Tyr1018 polyclonal antibody (Platelet-Derived Growth Factor Receptor Alpha) binds to endogenous Platelet-Derived Growth Factor Receptor Alpha only when phosphorylated at Tyr1018.
Immunogen Sequence	