

Anti-Phospho-PRKCD-Ser359 antibody (STJ196383)

STJ196383

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Protein Kinase C Delta Type-Ser359 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5580
Gene Symbol	PRKCD
Uniprot ID	KPCD_HUMAN
Immunogen	Synthesized phospho peptide around human PKC Delta (Ser359)
Immunogen Region	
Specificity	Phospho-PRKCD-Ser359 polyclonal antibody (Protein Kinase C Delta Type) binds to endogenous Protein Kinase C Delta Type only when phosphorylated at Ser359.
Immunogen Sequence	