

Anti-Phospho-PPIG-Ser376 antibody (STJ196389)

STJ196389

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Peptidyl-Prolyl Cis-Trans Isomerase G-Ser376 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9360
Gene Symbol	PPIG
Uniprot ID	PPIG_HUMAN
Immunogen	Synthesized phospho peptide around human PPIG (Ser376)
Immunogen Region	
Specificity	Phospho-PPIG-Ser376 polyclonal antibody (Peptidyl-Prolyl Cis-Trans Isomerase G) binds to endogenous Peptidyl-Prolyl Cis-Trans Isomerase G only when phosphorylated at Ser376.
Immunogen Sequence	