

Anti-Phospho-RAD50-Ser635 antibody (STJ196395)

STJ196395

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Dna Repair Protein Rad50-Ser635 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10111
Gene Symbol	RAD50
Uniprot ID	RAD50_HUMAN
Immunogen	Synthesized phosho peptide around human Rad50 (Ser635)
Immunogen Region	
Specificity	Phospho-RAD50-Ser635 polyclonal antibody (Dna Repair Protein Rad50) binds to endogenous Dna Repair Protein Rad50 only when phosphorylated at Ser635.
Immunogen Sequence	