

Anti-Phospho-RANBP3-Ser58 antibody (STJ196396)

STJ196396

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Ran-Binding Protein 3-Ser58 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	8498
Gene Symbol	RANBP3
Uniprot ID	RANB3_HUMAN
Immunogen	Synthesized phosho peptide around human RanBP3 (Ser58)
Immunogen Region	
Specificity	Phospho-RANBP3-Ser58 polyclonal antibody (Ran-Binding Protein 3) binds to endogenous Ran-Binding Protein 3 only when phosphorylated at Ser58.
Immunogen Sequence	