

Anti-Phospho-POLR2A-Thr4 antibody (STJ196408)

STJ196408

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Dna-Directed Rna Polymerase Ii Subunit Rpb1-Thr4 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5430
Gene Symbol	POLR2A
Uniprot ID	RPB1_HUMAN
Immunogen	Synthesized phosho peptide around human Rpb1 CTD (Thr4)
Region	
Specificity	Phospho-POLR2A-Thr4 polyclonal antibody (Dna-Directed Rna Polymerase Ii Subunit Rpb1) binds to endogenous Dna-Directed Rna Polymerase Ii Subunit Rpb1 only when phosphorylated at Thr4.
Immunogen Sequence	