

Anti-Phospho-SCRIB-Ser1220 antibody (STJ196413)

STJ196413

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Protein Scribble Homolog-Ser1220 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23513
Gene Symbol	SCRIB
Uniprot ID	SCRIB_HUMAN
Immunogen	Synthesized phosho peptide around human Scribble (Ser1220)
Immunogen Region	
Specificity	Phospho-SCRIB-Ser1220 polyclonal antibody (Protein Scribble Homolog) binds to endogenous Protein Scribble Homolog only when phosphorylated at Ser1220.
Immunogen Sequence	