

Anti-Phospho-SEMA4B-Ser825 antibody (STJ196414)

STJ196414

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Semaphorin-4B-Ser825 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10509
Gene Symbol	SEMA4B
Uniprot ID	SEM4B_HUMAN
Immunogen	Synthesized phospho peptide around human Semaphorin 4B (Ser825)
Immunogen Region	
Specificity	Phospho-SEMA4B-Ser825 polyclonal antibody (Semaphorin-4B) binds to endogenous Semaphorin-4B only when phosphorylated at Ser825.
Immunogen Sequence	