

Anti-Phospho-SHC1-Tyr317 antibody (STJ196418)

STJ196418

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Shc-Transforming Protein 1-Tyr317 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6464
Gene Symbol	SHC1
Uniprot ID	SHC1_HUMAN
Immunogen	Synthesized phosho peptide around human Shc (Tyr317)
Immunogen Region	
Specificity	Phospho-SHC1-Tyr317 polyclonal antibody (Shc-Transforming Protein 1) binds to endogenous Shc-Transforming Protein 1 only when phosphorylated at Tyr317.
Immunogen Sequence	