

## Anti-Phospho-INPPL1-Tyr986/987 antibody (STJ196420)

STJ196420

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Phospho-Phosphatidylinositol 3-4-5-Trisphosphate 5-Phosphatase 2-Tyr986/987 is suitable for use in Western Blot research applications.
<b>Applications</b>	WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Rat, Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
<b>Dilution Range</b>	WB 1:1000-2000
<b>Formulation</b>	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">3636</a>
<b>Gene Symbol</b>	<a href="#">INPPL1</a>
<b>Uniprot ID</b>	<a href="#">SHIP2_HUMAN</a>
<b>Immunogen</b>	Synthesized phosho peptide around human SHIP2 (Tyr986 and 987)
<b>Region</b>	
<b>Specificity</b>	Phospho-INPPL1-Tyr986/987 polyclonal antibody (Phosphatidylinositol 3-4-5-Trisphosphate 5-Phosphatase 2) binds to endogenous Phosphatidylinositol 3-4-5-Trisphosphate 5-Phosphatase 2 only when phosphorylated at Tyr986/987.
<b>Immunogen Sequence</b>	