

Anti-Phospho-SQSTM1-Ser349 antibody (STJ196432)

STJ196432

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Sequestosome-1-Ser349 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8878
Gene Symbol	SQSTM1
Uniprot ID	SQSTM_HUMAN
Immunogen	Synthesized phosho peptide around human SQSTM1 (Ser349)
Immunogen Region	
Specificity	Phospho-SQSTM1-Ser349 polyclonal antibody (Sequestosome-1) binds to endogenous Sequestosome-1 only when phosphorylated at Ser349.
Immunogen Sequence	