

Anti-Phospho-NCOA3-Thr24 antibody (STJ196437)

STJ196437

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Nuclear Receptor Coactivator 3-Thr24 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8202
Gene Symbol	NCOA3
Uniprot ID	NCOA3_HUMAN
Immunogen	Synthesized phospho peptide around human SRC-3 (Thr24)
Immunogen Region	
Specificity	Phospho-NCOA3-Thr24 polyclonal antibody (Nuclear Receptor Coactivator 3) binds to endogenous Nuclear Receptor Coactivator 3 only when phosphorylated at Thr24.
Immunogen Sequence	