

Anti-Phospho-STAT3-Tyr705 antibody (STJ196441)

STJ196441

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-Phospho-Signal Transducer And Activator Of Transcription 3-Tyr705 is suitable for use in Western Blot research applications. |
| Applications | WB |
| Host/Source | Rabbit |
| Reactivity | Human, Mouse, Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography. |
| Dilution Range | WB 1:1000-2000 |
| Formulation | PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide. |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 6774 |
| Gene Symbol | STAT3 |
| Uniprot ID | STAT3_HUMAN |
| Immunogen | Synthesized phosho peptide around human Stat3 (Tyr705) |
| Region | |
| Specificity | Phospho-STAT3-Tyr705 polyclonal antibody (Signal Transducer And Activator Of Transcription 3) binds to endogenous Signal Transducer And Activator Of Transcription 3 only when phosphorylated at Tyr705. |
| Immunogen Sequence | |