

Anti-Phospho-STAT5B-Tyr694 antibody (STJ196442)

STJ196442

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Signal Transducer And Activator Of Transcription 5b-Tyr694 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6777
Gene Symbol	STAT5B
Uniprot ID	STA5B_HUMAN
Immunogen	Synthesized phosho peptide around human Stat5 (Tyr694)
Region	
Specificity	Phospho-STAT5B-Tyr694 polyclonal antibody (Signal Transducer And Activator Of Transcription 5b) binds to endogenous Signal Transducer And Activator Of Transcription 5b only when phosphorylated at Tyr694.
Immunogen Sequence	