

Anti-Phospho-TOP2A-Ser1469 antibody (STJ196455)

STJ196455

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Dna Topoisomerase 2-Alpha-Ser1469 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:1000-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7153
Gene Symbol	TOP2A
Uniprot ID	TOP2A_HUMAN
Immunogen	Synthesized phosho peptide around human Topoisomerase II Alpha (Ser1469)
Immunogen Region	
Specificity	Phospho-TOP2A-Ser1469 polyclonal antibody (Dna Topoisomerase 2-Alpha) binds to endogenous Dna Topoisomerase 2-Alpha only when phosphorylated at Ser1469.
Immunogen Sequence	