

## Anti-Phospho-ULK1-Ser638 antibody (STJ196464)

STJ196464

### GENERAL INFORMATION

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Primary antibodies   |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-Phospho-Serine/Threonine-Protein Kinase Ulk1-Ser638 is suitable for use in Western Blot research applications. |
| <b>Applications</b>      | WB   |
| <b>Host/Source</b>       | Rabbit   |
| <b>Reactivity</b>        | Human, Mouse   |

### PRODUCT PROPERTIES

|                            |  |
|----------------------------|--|
| <b>Clonality</b>           | Polyclonal   |
| <b>Clone ID</b>            |  |
| <b>Concentration</b>       | 1 mg/mL  |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Purification</b>        | The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.          |
| <b>Dilution Range</b>      | WB 1:1000-2000   |
| <b>Formulation</b>         | PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.  |
| <b>Isotype</b>             | IgG  |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                           |   |
|---------------------------|---|
| <b>Gene ID</b>            | <a href="#">8408</a>  |
| <b>Gene Symbol</b>        | <a href="#">ULK1</a>  |
| <b>Uniprot ID</b>         | <a href="#">ULK1_HUMAN</a>  |
| <b>Immunogen</b>          | Synthesized phospho peptide around human ULK1 (Ser638)  |
| <b>Immunogen Region</b>   |   |
| <b>Specificity</b>        | Phospho-ULK1-Ser638 polyclonal antibody (Serine/Threonine-Protein Kinase Ulk1) binds to endogenous Serine/Threonine-Protein Kinase Ulk1 only when phosphorylated at Ser638. |
| <b>Immunogen Sequence</b> |   |

---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081