

## Anti-Phospho-ULK1-Ser757 antibody (STJ196465)

STJ196465

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Phospho-Serine/Threonine-Protein Kinase Ulk1-Ser757 is suitable for use in Western Blot research applications.
<b>Applications</b>	WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
<b>Dilution Range</b>	WB 1:1000-2000
<b>Formulation</b>	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	8408
<b>Gene Symbol</b>	ULK1
<b>Uniprot ID</b>	ULK1_HUMAN
<b>Immunogen</b>	Synthesized phospho peptide around human ULK1 (Ser757)
<b>Immunogen Region</b>	
<b>Specificity</b>	Phospho-ULK1-Ser757 polyclonal antibody (Serine/Threonine-Protein Kinase Ulk1) binds to endogenous Serine/Threonine-Protein Kinase Ulk1 only when phosphorylated at Ser757.
<b>Immunogen Sequence</b>	

---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081