

Anti-Thymidine Kinase antibody (STJ196479)

STJ196479

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Thymidine kinase is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Rabbit
Reactivity	HHV-1

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthesized peptide derived from human Thymidine kinase
Immunogen Region	
Specificity	Thymidine Kinase polyclonal antibody (Thymidine kinase) binds to endogenous Thymidine kinase.
Immunogen Sequence	