

Anti-Serine/threonine-protein kinase BGLF4 antibody (STJ196485)

STJ196485

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Serine/threonine-protein kinase BGLF4 is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Rabbit
Reactivity	HHV-4

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthesized peptide derived from human Serine/threonine-protein kinase BGLF4
Immunogen	
Region	
Specificity	Serine/threonine-protein kinase BGLF4 polyclonal antibody (Serine/threonine-protein kinase BGLF4) binds to endogenous Serine/threonine-protein kinase BGLF4.
Immunogen	
Sequence	