

Anti-Serine/threonine protein kinase UL97 antibody (STJ196487)

STJ196487

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Serine/threonine protein kinase UL97 is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Rabbit
Reactivity	HHV-5

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthesized peptide derived from human Serine/threonine protein kinase UL97
Immunogen Region	
Specificity	Serine/threonine protein kinase UL97 polyclonal antibody (Serine/threonine protein kinase UL97) binds to endogenous Serine/threonine protein kinase UL97.
Immunogen Sequence	