

Anti-POLD4 antibody (STJ196496)

STJ196496

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Dna Polymerase Delta Subunit 4 is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-200 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57804
Gene Symbol	POLD4
Uniprot ID	DPOD4_HUMAN
Immunogen	Synthesized peptide derived from human POLD4
Immunogen Region	
Specificity	POLD4 polyclonal antibody (Dna Polymerase Delta Subunit 4) binds to endogenous Dna Polymerase Delta Subunit 4.
Immunogen Sequence	