

Anti-PKIG antibody (STJ196500)

STJ196500

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Camp-Dependent Protein Kinase Inhibitor Gamma is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-200 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	11142
Gene Symbol	PKIG
Uniprot ID	IPKG_HUMAN
Immunogen	Synthesized peptide derived from human PKIG
Immunogen Region	
Specificity	PKIG polyclonal antibody (Camp-Dependent Protein Kinase Inhibitor Gamma) binds to endogenous Camp-Dependent Protein Kinase Inhibitor Gamma.
Immunogen Sequence	