

Anti-WT1-AS antibody (STJ196505)

STJ196505

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Putative Wilms Tumor Upstream Neighbor 1 Gene Protein is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-200 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	NA
Gene Symbol	WT1-AS
Uniprot ID	WIT1_HUMAN
Immunogen	Synthesized peptide derived from human WIT1
Immunogen Region	
Specificity	WT1-AS polyclonal antibody (Putative Wilms Tumor Upstream Neighbor 1 Gene Protein) binds to endogenous Putative Wilms Tumor Upstream Neighbor 1 Gene Protein.
Immunogen Sequence	