

## Anti-GNG4 antibody (STJ196522)

STJ196522

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Guanine Nucleotide-Binding Protein G (I/G (S/G (O Subunit Gamma-4 is suitable for use in Immunohistochemistry and ELISA research applications.
<b>Applications</b>	IHC, ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
<b>Dilution Range</b>	IHC-P 1:50-200 ELISA 1:5000-20000
<b>Formulation</b>	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">2786</a>
<b>Gene Symbol</b>	<a href="#">GNG4</a>
<b>Uniprot ID</b>	<a href="#">GBG4_HUMAN</a>
<b>Immunogen</b>	Synthesized peptide derived from human GNG4
<b>Immunogen Region</b>	
<b>Specificity</b>	GNG4 polyclonal antibody (Guanine Nucleotide-Binding Protein G (I/G (S/G (O Subunit Gamma-4) binds to endogenous Guanine Nucleotide-Binding Protein G (I/G (S/G (O Subunit Gamma-4.
<b>Immunogen Sequence</b>	