

Anti-DEFB1 antibody (STJ196527)

STJ196527

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Beta-Defensin 1 is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-200 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1672
Gene Symbol	DEFB1
Uniprot ID	DEFB1_HUMAN
Immunogen	Synthesized peptide derived from human Beta-defensin 1
Immunogen Region	
Specificity	DEFB1 polyclonal antibody (Beta-Defensin 1) binds to endogenous Beta-Defensin 1.
Immunogen Sequence	