

## Anti-SCGB1D4 antibody (STJ196570)

STJ196570

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Secretoglobin Family 1d Member 4 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
<b>Applications</b>	IHC-P, IF-P
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Rat, Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
<b>Dilution Range</b>	IHC-P 1:50-200
<b>Formulation</b>	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">404552</a>
<b>Gene Symbol</b>	<a href="#">SCGB1D4</a>
<b>Uniprot ID</b>	<a href="#">SG1D4_HUMAN</a>
<b>Immunogen</b>	Synthesized peptide derived from human SG1D4
<b>Immunogen Region</b>	
<b>Specificity</b>	SCGB1D4 polyclonal antibody (Secretoglobin Family 1d Member 4) binds to endogenous Secretoglobin Family 1d Member 4.
<b>Immunogen Sequence</b>	