

Anti-PATE4 antibody (STJ196575)

STJ196575

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Prostate And Testis Expressed Protein 4 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-P, IF-P
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-200
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	399968
Gene Symbol	PATE4
Uniprot ID	PATE4_HUMAN
Immunogen	Synthesized peptide derived from human PATE4
Immunogen Region	
Specificity	PATE4 polyclonal antibody (Prostate And Testis Expressed Protein 4) binds to endogenous Prostate And Testis Expressed Protein 4.
Immunogen Sequence	