

Anti-SDHAF4 antibody (STJ196580)

STJ196580

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Succinate Dehydrogenase Assembly Factor 4-Mitochondrial is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-P, IF-P
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-200
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	135154
Gene Symbol	SDHAF4
Uniprot ID	SDHF4_HUMAN
Immunogen	Synthesized peptide derived from human CF057
Immunogen Region	
Specificity	SDHAF4 polyclonal antibody (Succinate Dehydrogenase Assembly Factor 4-Mitochondrial) binds to endogenous Succinate Dehydrogenase Assembly Factor 4-Mitochondrial.
Immunogen Sequence	