

Anti-TMSB15B antibody (STJ196586)

STJ196586

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Thymosin Beta-15b is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-P, IF-P
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-200
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11013/286527
Gene Symbol	TMSB15B
Uniprot ID	TB15B_HUMAN
Immunogen	Synthesized peptide derived from human TB15B
Immunogen Region	
Specificity	TMSB15B polyclonal antibody (Thymosin Beta-15b) binds to endogenous Thymosin Beta-15b.
Immunogen Sequence	