

Anti-FBXL5 antibody (STJ196587)

STJ196587

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-F-Box/Lrr-Repeat Protein 5 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	26234
Gene Symbol	FBXL5
Uniprot ID	FBXL5_HUMAN
Immunogen	Synthesized peptide derived from human FBXL5
Immunogen Region	
Specificity	FBXL5 polyclonal antibody (F-Box/Lrr-Repeat Protein 5) binds to endogenous F-Box/Lrr-Repeat Protein 5.
Immunogen Sequence	