

Anti-LINC02877 antibody (STJ196592)

STJ196592

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Putative Uncharacterized Protein Encoded By Linc02877 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-P, IF-P
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-200
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	NA
Gene Symbol	LINC02877
Uniprot ID	CC079_HUMAN
Immunogen	Synthesized peptide derived from human CC079
Immunogen Region	
Specificity	LINC02877 polyclonal antibody (Putative Uncharacterized Protein Encoded By Linc02877) binds to endogenous Putative Uncharacterized Protein Encoded By Linc02877.
Immunogen Sequence	